

Generation of a Standard Curve to Convert CTG Values to Cell Number for PC9 Cells

- Protocol
 - Seed PC9 cells at different cell numbers per well in triplicate in 90ul RPMI+10%FBS growth medium
 - Incubate at 37C, 5% CO2
 - 24h after seeding, image wells by IncuCyte, then read wells by CTG
- Analysis
 - For all CTG values from wells containing cells, subtract CTG value from empty well (raw data minus background)
 - Calculate average value and standard deviation for each triplicate
 - Plot cell number seeded (X axis) vs. CTG reading (Y axis)

Plate Map - 1000 to 64000 cells/well by 2X increments, done in triplicate

[illegible]

Day 1 IncuCyte Image

Cell \ well

1000

2000








4000

8000

16000


32000

64000




RLU =

- CTG reading from wells with cells
- Subtract empty well CTG value
- Average of triplicates



- Use standard curve to convert CTG values to cell numbers from all subsequent drug treatment experiments (formula at CTG cell# on excel)
- Use 2000 cells as the starting number of cells at the initiation of drug treatment (day 0 baseline)

Formula at CTG cell# on excel



Single Agent and Combo Effects of Dacomitinib and Osimertinib on PC9 Cell Viability in 3 Day Assays

- Protocol
 - Seed PC9 cells at 2000 cells per well in 80 or 90ul RPMI+10%FBS growth medium
 - Incubate at 37C, 5% CO2
 - 24h after seeding (day 0 baseline), add single agent or combo drug treatment as indicated in plate map, with duplicate wells for each condition
 - After 3 days of drug treatment (day 3), read wells by CTG
- Bar Graph Analysis
 - For all CTG values from wells containing cells, subtract CTG value from empty well (raw data minus background)
 - Convert all CTG values to cell numbers using standard curve (formula at CTG cell# on excel)
 - Bar plot of cell numbers per selected drug treatment conditions in duplicate
 - Use 2000 cells seeded CTG read at day 0 converted through formula to cell number as untreated starting cell number (day 0)
- Concentration – Response Curve Analysis
 - All data is from CTG day 3 read, plates 3 and 4, drug treatment and control DMSO vehicle
 - For all CTG values, subtract empty well background
 - $\% \text{ control} = 100\% \times (\text{treatment CTG read} / \text{DMSO CTG read})$
 - Concentration-response curve plotted log drug concentration (X-axis) versus % control (Y-axis) using GraphPad PRISM



Plate map



Plate 3: seed 2K cells/well; start drug treatment next day; CTG read after 3 days drug treatment

[illegible]




Plate 4: seed 2K cells/well; start drug treatment next day; CTG read after 3 days drug treatment

[illegible]

Bar Graph Plots



Concentration-Response Curves



PF804=Dacomitinib
PF9696=Osimertinib